

iSupport Platform

For your daily care

Team members:

Burdisso Sara s275424

Centonze Ilenia s277389

Ferrara Mauro s276932


Pasotti Andrea s275100



The idea - overview

Objective:

To help clinicians and caregivers follow individuals in depressive states during their treatment period by providing a remote monitoring and support platform.



Protocolli di Comunicazione

- *REST Web Service*

It allows to synchronize the communication between the parties.

Catalog: accesses the system information present in a JSON file and communicates it to the other actors;


ThingSpeak: database with data coming from the sensors. The data is accessed with a GET request and provided via a message in JSON format.

- *MQTT*


Lightweight and flexible protocol suitable for IoT applications.

- **Publish** messages in JSON format using QoS 2
- **Receive** messages using QoS 2
- **Message topic:** hierarchical, containing patient information

Use Case Diagram - proposal



Use Case Diagram - modified



Catalog

- It communicates with all the other actors in the platform exploiting REST communication.



```
{  
  "Catalog_url":"http://127.0.0.1:8080",  
  "baseTopic":"iSupport/",  
  "LightShift":{  
    "clientID" : "LigthShiftMS",  
    "endTopic":["actuators/Light"]  
  },  
  ...  
}
```

Catalog


- Catalog.json contains information about general configuration, registered devices and registered patients.



```
{
  "broker": {
    "IPAddress": "test.mosquitto.org",
    "port": 1883
  },
  "token": "1711013176:AAH6bCdVA0d0dL8fX4FIfCPGEqpEFNjN1d4",
  "deviceList": [
    {
      "deviceName": "smartwatch_1",
      "deviceID": "1",
      "measureType": [
        "HeartRate",
        "Accelerometer"
      ],
      "availableServices": [
        "MQTT"
      ],
      "servicesDetails": [
        {
          "serviceType": "MQTT",
          "topic": [
            "iSupport/1/sensors/Body"
          ]
        }
      ],
      "lastUpdate": 1621764472.704841
    }
  ],
  "patientList": [
    {
      "patientID": "1",
      "patientName": "Mario",
      "patientSurname": "Rossi",
      "LightsSchedule": "7:30-8:30",
      "connectedDevices": [
        {
          "measure": [
            "Heart Rate",
            "Acceration"
          ],
          "deviceID": "smartwatch_1"
        },
        {
          "telegramIDs": [1234567],
          "thingspeakInfo": {
            "apikeys": [
              "ZQ2AJZABDSC4MXXG",
              "SV1HOET58CE7225M"
            ]
          },
          "channel": 1342741
        }
      ],
      "lastUpdate": 1622981412.5148401
    }
  ]
}
```


Catalog

- It communicates with all the other actors in the platform exploiting REST communication.
- It has 3 main functions:
 - Retrieve (**GET**)
 - Add (POST)
 - Update (PUT)




Catalog

- It communicates with all the other actors in the platform exploiting REST communication.
- It has 3 main functions :
 - Retrieve (GET)
 - Add (**POST**)
 - Update (PUT)



Catalog

- It communicates with all the other actors in the platform exploiting REST communication.
- It has 3 main functions:
 - Retrieve (GET)
 - Add (POST)
 - Update (**PUT**)




```
Body={
  "patientID": «1»,
  "patientName": "Francesco",
  "patientSurname": "Bianchi",
  "LightsSchedule": "17:00-18:30",
  "connectedDevices": [
    {
      "measure": "Temperature",
      "deviceID": "t_1"
    }
  ],
  "telegramIDs": [],
  "thingspeakInfo": {
    "apikey": [
      "ZQ2AJZABDSC4MXXG",
      "SV1HOET58CE7225M"
    ],
    "channel": 1342741
  },
  "lastUpdate": ""
}
```

```
Body={
  "deviceName": "smartwatch_1",
  "deviceID": "1",
  "measureType": [
    "HeartRate",
    "Accelerometer"
  ],
  "availableServices": [
    "MQTT"
  ],
  "servicesDetails": [
    {
      "serviceType": "MQTT",
      "topic": [
        "iSupport/1/sensors/Body"
      ]
    }
  ],
  "lastUpdate": 1623097970.3872042
}
```

Catalog

- Catalog keeps updated the "deviceList" with the devices correctly connected and removes all the devices that are disconnected.

```
"deviceList": [  
  {  
    "deviceName": "smartwatch_1",  
    "deviceID": "1",  
    "measureType": [  
      "HeartRate",  
      "Accelerometer"  
    ],  
    "availableServices": [  
      "MQTT"  
    ],  
    "servicesDetails": [  
      {  
        "serviceType": "MQTT",  
        "topic": [  
          "iSupport/1/sensors/Body"  
        ]  
      }  
    ],  
    "lastUpdate": 1621764472.704841  
  }  
]
```



Device Connector

Devices implemented in the simulation:

- **Temperature and humidity sensors**
- **Motion Sensor**
- Heart rate sensor
- Accelerometer sensor
- Actuator for air Conditionair system
- Smart light bulb actuator

- Temperature and Humidity:

Extraction of a value from uniform distributions that have specific limits that change according to the month in which publication occurs. User sets mode of generation:

- in the suitable range
- out of range

- Motion Sensor:

It is a value (1/0) set by the user

Device Connector


Devices implemented in the simulation:

- Temperature and humidity sensors
- Motion Sensor
- **Heart rate sensor**
- **Accelerometer sensor**
- Actuator for air Conditionair system
- Smart light bulb actuator

-HR and Accelerometer measurments:

Accelerometer measurments are obtained from 2 database:
REST.txt and *SPORT.txt*

• REST configuration



+

REST.txt

0.06275
0.04957
0.11802
0.28902
0.10596
0.556840...

Device Connector


Devices implemented in the simulation:

- Temperature and humidity sensors
- Motion Sensor
- **Heart rate sensor**
- **Accelerometer sensor**
- Actuator for air Conditionair system
- Smart light bulb actuator

-HR and Accelerometer measurments:

Accelerometer measurments are obtained from 2 database:
REST.txt and *SPORT.txt*

• Sport configuration



+

SPORT.txt

6.920512
7.049275
7.176033
7.291157
7.398414
8.597905
8.611896...

Device Connector


Devices implemented in the simulation:

- Temperature and humidity sensors
- Motion Sensor
- **Heart rate sensor**
- **Accelerometer sensor**
- Actuator for air Conditionair system
- Smart light bulb actuator

-HR and Accelerometer measurments:

Accelerometer measurments are obtained from 2 database:
REST.txt and *SPORT.txt*

- «Danger» configuration



+

REST.txt

0.06275

0.04957

0.11802

0.28902


0.10596

0.556840...



Device Connector

Devices implemented in the simulation:

- Temperature and humidity sensors
 - Presence Sensor
 - Heart rate sensor
 - Accelerometer sensor
 - **Actuator for air Conditionair system**
 - **Smart light bulb actuator**
- 

Actuation Commands:

- Switching on the air Conditionair system to maintain optimal temperature and humidity levels.
- Switching on smart lights to enable luminothrapy.

Device Connector

Communication with the other actors is done using **SenML format**.

```
self.__message={
    "bn": 1,
    "e": [{
        "n": "Heart Rate",
        "value": 61,
        "timestamp": "123456",
        "u": "bpm"
    },
    {
        "n": "Accelerometer",
        "value": 0.060369,
        "timestamp": "123456",
        "u": "m/s2"
    }
    ]
}
```

Communication with

- **Catalog**: to register devices or refresh registration with a REST communication.
- **ThingSpeak Adaptor**: which receives all measurements from the sensors published by Device Connector.
- **Home system Control**: device connector publishes sensor measurements to this microservice, which sends actuation commands to switch on and off the temperature control system.
- **Light Shift**: from which device connector receives actuation command to switch on/off smart light bulbs for luminotherapy.

Patient Control



Goal: Monitoring the patient and detecting of panic attack

How it works

- Analyze heart rate and acceleration data

Communication with


- ThingSpeak: for getting measurement of 10 minutes windows.



Patient Control

How it works

- Analyze heart rate and acceleration data



- Communicate an alert in case of panic attack to the doctor and caregiver

Communication with


- ThingSpeak: for getting measurement of 10 minutes windows.
- ThingSpeak: for recording the panic attack event
- Telegram bot: for communicating the panic attack event

Home system Control

Goal: Monitoring home for having a comfortable enviroment

How it works


- Manages temperature and humidity inside the house
- Verifies patient presence in the home
- Avoid wasting energy by establishing activation time slots



Home system Control

How it works

- Manages temperature and humidity inside the house
- Verifies patient presence in the home
- Avoid wasting energy by establishing activation time slots



Home system Control

How it works

- Manages temperature and humidity inside the house
- Verifies patient presence in the home
- Avoid wasting energy by establishing activation time slots

Communication with

- Device connector:
 - Continuously receives data from Temperature and Humidity sensors
 - Sends actuation commands if it is identified a condition to switch ON or OFF the Air Conditioning System

Light Shift



How it works

- Manages the lighting system in relation to the physician's established schedule for luminotherapy

Communication with


- Catalog:
 - to obtain the times in which to activate the luminotherapy
- Device Connectors:
 - to send actuations commands to turn the lights ON or OFF








ThingSpeak

It's a platform used like a database in which there are all measurement got from different sensors implemented on iSupport platform

It uses an adaptor *ThingSpeak Adaptor* which:

- Allows you to implement communication via MQTT.
 - Increases the scalability of the platform
 - Increases interoperability between components.
- 

ThingSpeak third part Application User Interface



ThingSpeak

It's a platform used like a database in which there are all measurement got from different sensors implemented on iSupport platform

It uses an adaptor *ThingSpeak Adaptor* which:

- Allows you to implement communication via MQTT.
- Increases the scalability of the platform
- Increases interoperability between components.

Communication with

- Device Connector: continuous updating of measurements
- Data Analysis: sending measurements received from the device connector.
- Patient Control: sending measurements related to specific observation windows.


Data Analysis

How it works

- Process patient measurements to get weekly or monthly statistics:
 - **Weekly report about average heart rate in 3 time bands.**
 - Weekly report about patient's activity status (percentage of low, medium, high intensity obtained as comparing rate and acceleration measurements)
 - Number of monthly panic attacks
 - Weekly report about how much time patient spends in the bedroom (level of sedentary lifestyle)

Communication with

- ThingSpeak Adaptor: to get measurements
- Node-RED: to create statistical graphs



Node-RED interface


Data Analysis

How it works

- Process patient measurements to get weekly or monthly statistics:
 - Weekly report about average frequency values in 3 time bands.
 - **Weekly report about patient's activity status (percentage of low, medium, high intensity obtained as comparing rate and acceleration measurements.**
- Number of monthly panic attacks.
- Weekly report about how much time patient spends in the bedroom (level of sedentary lifestyle).

Communication with

- ThingSpeak Adaptor: to get measurements
- Node-RED: to create statistical graphs



Node-RED interface


Data Analysis

How it works

- Process patient measurements to get weekly or monthly statistics:
 - Weekly report about average frequency values in 3 time bands.
 - Weekly report about patient's activity status (percentage of low, medium, high intensity obtained as comparing rate and acceleration measurements).
 - **Number of monthly panic attacks.**
 - **Weekly report about how much time patient spends in the bedroom (level of sedentary lifestyle).**

Communication with

- ThingSpeak Adaptor: to get measurements
- Node-RED: to create statistical graphs




Node-RED interface

Node-RED

How it works

- Nodes:


- MQTT communication node
- *Python function*
- Two blocks to visualize msg.payload in debug *display*
- A *rbe* block
- A *digital display* and three *chart*





Telegram Bot


- Emergency management
- Allow the doctor and the caregiver to always be updated on the Heart Rate, Movement and Recurrence trends of panic attacks.

- '/start' function:
 - Telegram user registration on platform
 - '/status' function:
 - Visualization of ThingSpeak graphs about vital signs
 - Emergency management:
 - Getting alert message from the Patient Control microservice in case of Panic Attack.
- 

Telegram Bot


- Emergency management
- Allow the doctor and the caregiver to always be updated on the Heart Rate, Movement and Recurrence trends of panic attacks.

- '/start' function:
 - Telegram user registration on platform




Telegram Bot

- Emergency management
- Allow the doctor and the caregiver to always be updated on the Heart Rate, Movement and Recurrence trends of panic attacks.



- '/status' function:


- Visualization of ThingSpeak graphs about vital signs



Telegram Bot


- Emergency management
- Allow the doctor and the caregiver to always be updated on the Heart Rate, Movement and Recurrence trends of panic attacks.

- Emergency management:
 - Getting alert message from the Patient Control microservice in case of Panic Attack.





Telegram Bot

- Emergency management
 - Allow the doctor and the caregiver to always be updated on the Heart Rate, Movement and Recurrence trends of panic attacks.
- 

Communication with

- Patient Control: which publishes a message to the telegram bot if it detects a panic attack
- ThingSpeak: required for viewing graphs
- Catalog: retrieve chatId and token



Thanks for your attention

Team 3:

Burdisso Sara s275424

Centonze Ilenia s277389

Ferrara Mauro s276932

Pasotti Andrea s275100